



Removable part for
massage position.
جزء قابل للفك لوضع المساج.

BMI
BORG EL ARAB
MEDICAL INDUSTRIES



Sand blasting
تنظيف الحديد بالترميز



Co2 Welding
لحام ثاني اكسيد الكربون

Optional

LINAK

colors



Electric examination bed

- Luxurious leather mattress with high density foam, divided to 2 parts (backrest, seat) removable part for massage position.
- Leather color can be changed per request.
- Bed frame is made of 3 x 5 steel sections, 1.8mm thickness, electrostatic coated.
- Co2 Welding.
- Sand blasting.
- Dimensions: Length: 195cm Width: 60cm Height adjustment from 50-75cm.
- Positions: Back rest lifting from 0° to 70°.
- 2 motors, noiseless.
- Safe working load (SWL) 200kg.
- Optional: LINAK motors (Denmark).

EXM 2E30

سرير كشف كهربائي

- مرتبة فاخرة من الجلد عالية الكثافة مقسمة إلي جزئين (مسند الظهر , المقعد) جزء قابل للفك لوضع المساج .
- إمكانية تغير لون الجلد حسب الطلب .
- جسم السرير مصنوع من علب 3 × 5 سم السمك 1.8 مم . دهان إلكتروستاتك .
- لحام ثاني أكسيد الكربون .
- تنظيف الحديد بالترميز .
- الأبعاد : الطول 195سم العرض : 60 سم .
- الارتفاع من 50سم إلي 75سم .
- الحركات : زاوية حركة الظهر من 0° إلي 70° .
- 2 موتور , صوت منخفض .
- وحدة تحكم لليد .
- الحمل الأمان للتشغيل 200 كيلوجرام .
- اختياري : إمكانية تركيب مولدات ماركه LINAK (صناعة دنيماركي) .