

ET-85001A



Emergency trolley

- ABS body, metal ball - bearing rails for drawers, Aluminum columns
- 2 small & 3 big drawers with inner partitions can be organized
- Accessories: needle disposal holder two dust baskets, power outlet 5m wire CPR board, oxygen tank holder Monitor shelf, IV pole, files holder, Medicin preparation rack.
- 5 inches castors two of them with brakes
- Centralized lock for drawers.
- Stainless steel guard rail.
- Concave top cover transparent. soft plastic on top board
- Size: 850*520*950mm
- Optional: Size: 750*475*930 mm with 4 inches castors

تولاي طوارئ

- الجسم مصنوع من ABS، مجاري الأدرج مزودة برومان بلي، القوائم مصنوعة من الألومنيوم.
- 2 درج صغير و3 درج كبير، فواصل متحركة من الداخل يمكن تنظيمها.
- الأكسسورات: سلة للتخلص من الإبر، 2 سلة للنفايات
- مخرج كهرباء (سلك 5 مم)، لوحة إنعاش قلبي، حامل إسطوانة أكسجين، رف لجهاز المونيتور، حامل محلول، حامل للملفات، رف تحضير الأدوية.
- عجل 5 بوصة، 2 منهم بفرامل.
- قفل مركزي للأدرج.
- حاجز حماية من الستائلس.
- غطاء بلاستيكي شفاف علي السطح العلوي.
- الأبعاد: (850*520*950مم).
- إختياري: الأبعاد: (930*475*750مم) بعجل 4 بوصة