

colors



Pediatric examination bed

EXB 1M50-1

سرير كشف أطفال

- Luxurious leather mattress with high density foam, divided to 2 parts (backrest, seat).
- Leather color can be changed per request.
- Any photo or drawing can be printed on the bed sides.
- Bed frame is made of 3 x 3 steel sections, 1.8mm thickness, electrostatic coated bed sides made of thermosetting laminated plastic sheet (HPL).
- Dimensions: Length: 150 cm, Width: 60 cm, Height 80 cm
- Positions: Back rest lifting from 0° to 70°.
- Storage space under the bed.
- Safe working load (SWL) 100kg.

- مرتبة فاخرة من الجلد عالية الكثافة . مقسمة إلى جزئين (مسند الظهر , المقعد) . إمكانية تغير لون الجلد حسب الطلب .
- إمكانية طباعة أي صورة أو رسمة على جوانب السرير .
- جسم السرير مصنوع من علب حديد 3 × 3 سم السمك 1.8 مم ، دهان إلكتروستاتك .
- جوانب السرير مصنوع من شرائح البلاستيك المضغوطة (HPL) .
- الأبعاد : الطول 150 سم ، العرض : 60 سم ، الارتفاع 80 سم .
- الحركات : زاوية حركة الظهر من 0° إلى 70° .
- مساحة للتخزين أسفل السرير .
- الحمل الأمان للتشغيل 100 كيلوجرام .