



Emergency trolley

- HPL body, metal ball - bearing rails for drawers, Aluminum columns
- 3 small & 2 big drawers with inner partitions can be organized
- Accessories: needle disposal holder two dust baskets, power outlet 5m wire CPR board, oxygen tank holder Monitor shelf, IV pole, files holder, Medicin preparation rack.
- 3 inches castors with brakes .
- Centralized lock for drawers.
- Stainless steel guard rail.
- Concave top cover transparent. soft plastic on top board
- Size: 850*520*950mm

Model F-CRS-0014

ترولي طوارئ

- الجسم مصنوع من HPL, مجاري الأدرج مزودة برومان بلي, القوائم مصنوعة من الألومنيوم.
- 3 درج صغير و2 درج كبير, فواصل متحركة من الداخل يمكن تنظيمها.
- الإكسسورات: سلة للتخلص من الإبر, 1 سلة للنفايات
- مخرج كهرباء (سلك 5 مم), لوحة إنعاش قلبي
- حامل إسطوانة أكسجين, رف لجهاز المونيتور
- حامل محلول, حامل للملفات, رف تحضير الأدوية
- 4 عجل 3 بوصة بفرامل.
- قفل مركزي للأدرج.
- حاجز حماية من الستانلس.
- غطاء بلاستيك شفاف علي السطح العلوي.
- الأبعاد: (950*520*850مم).